



LX Compression System

High-quality, professional compression machines

The Doctor Life massage systems are developed and manufactured by Daesung Maref Ltd. This company was founded in 1986 in South Korea and has become one of the largest suppliers of pneumatic massage systems for medical, wellness and sports applications.

The company has a Research & Development department, where developers work with clinicians to create the best possible equipment. The systems are exported to over 70 countries and have been awarded the European certificate for medical devices CE0499 and the US FDA approval certificate.



*Visit meridiusmedical.com or call 01392 829133
for more information on Meridius products*




LX7 & LX9 Lymphedema Therapy

LX Pump Specifications:

- Adjustable Pressure Range: 20 -200 mmHg
- Safe Pressure Setting LOCK Function
- Illuminated Start/Pause/Stop Button
- Cycle Time: 60 seconds
- Therapy Timer: 10, 20, 30 minutes
- Chamber Control Options: On/FOCUS, Off/SKIP
- Digital Interface, Controller & Pressure Sensors
- DME Warranty: 2 yrs. (pump), 6 mon. (sleeves)
- Device Size: 10" x 5" x 7" (2.2 kg)

LX Sleeve Types:


- Legs: Large - 28"; XL - 32"; XXL - 38"
- Arms
- Waist
- Center Body
- Expanders for added width
- Integrated Hose
- Convenient One-Touch Connectors



Treatment Mode A - Peristaltic

Peristaltic treatment could help with:

- Mixed Arterial-Venous Disease with Edema
- DVT Prophylaxis
- Sensitive areas in the foot or ankle



Treatment Mode B - Sequential

Sequential treatment could help with:

- Lymphedema
- Venous Disease with Edema and/or wounds
- Inflammation Control
- Muscle Recovery