

PRODUCT CATALOGUE

2021 / 2022



Meridius Medical is based in Oberhausen, Germany. We are the sister company of AcuPrime, which is based in Exeter, United Kingdom.

Since 2001, we have been aspiring to grow our company by providing our services to the professional community.

With the meridius range, we strive to improve interconnectivity between the world and the body. We want to encourage people to live well.

We are the exclusive European distributor of Dongbang as well as the exclusive UK distributor of Spoband and Dr. Life intermittent pneumatic compression devices.

Our aim is to bring "Excellence at Sensible Prices" to healthcare professionals, helping restore vitality and wellbeing in both individuals and society.

www.eu.meridiusmedical.com
eu@meridiusmedical.com
+49 (0) 208 20580580

01

Acupuncture Needles

18

Physiotherapy

25

Cupping Therapy & Massage

36

Press Pellet & Probes

46

Electrical Devices

50

Moxa & Heating Devices

68

Treatment Room Equipment

meridius® Acupuncture Needles

meridius® Jet



Premium Plastic Handled Needles with Guide Tube

- One needle with a tube per blister
- Premium plastic handle for better grip
- Packed in 100 pcs per box
- EO gas sterilized

Code	Size
MEJET1215	0.12 x 15 mm
MEJET1230	0.12 x 30 mm
MEJET1415	0.14 x 15 mm
MEJET1430	0.14 x 30 mm
MEJET1615	0.16 x 15 mm
MEJET1630	0.16 x 30 mm
MEJET1815	0.18 x 15 mm
MEJET1830	0.18 x 30 mm
MEJET2015	0.20 x 15 mm
MEJET2030	0.20 x 30 mm
MEJET2530	0.25 x 30 mm
MEJET2540	0.25 x 40 mm
MEJET3030	0.30 x 30 mm
MEJET3040	0.30 x 40 mm
MEJET3050	0.30 x 50 mm

meridius® Eco



Stainless Steel Needles with Korean-style Handles

- Ten needles with a tube per pouch
- The needle features a Korean-style handle
- Packed in 50 pouches (500 needles) per box
- Gamma ray sterilized for a 5-year shelf life

Code	Size
MECO2030	0.20 x 30 mm
MECO2040	0.20 x 40 mm
MECO2530	0.25 x 30 mm
MECO2540	0.25 x 40 mm
MECO3030	0.30 x 30 mm
MECO3040	0.30 x 40 mm
MECO3050	0.30 x 50 mm

Acupuncture Needles

meridius® Acupuncture Needles

meridius® Jet



Premium Plastic Handled Needles without a Guide Tube

- One needle per blister
- Premium plastic handle for better grip
- Packed in 100 pcs per box
- EO gas sterilized

[illegible]

meridius® Eco



Stainless Steel Needles with Plastic Handles

- Five needles and one guide tube per blister
- The needle features a premium plastic handle
- 500 needles per box
- Sterilised using EO gas, with a 5 year shelf life.

Code	Size
OA5GME1615	0.16 x 15mm
OA5GME1630	0.16 x 30mm
OA5GME2030	0.20 x 30 mm
OA5GME2530	0.25 x 30 mm
OA5GME2540	0.25 x 40mm
OA5GME3050	0.30 x 50 mm

Acupuncture Needles

meridius® Pure



Silicone-Free Needles with Copper Handles

- One silicone-free needle without a tube per blister
- Copper handle for electrical and thermal conduction
- Packed in 100 pcs per box
- EO gas sterilized

Code	Size
MEPU2015	0.20 x 15 mm
MEPU2025	0.20 x 25 mm
MEPU2525	0.25 x 25 mm
MEPU2530	0.25 x 30 mm
MEPU3030	0.30 x 30 mm
MEPU3540	0.35 x 40 mm
MEPU4030	0.40 x 30 mm

meridius® Renew



Highly Visible Plastic Handled Needles without Tubes

- Five needles without tubes per blister
- Colour-coded plastic handles for better visibility
- Packed in 100 pcs per box
- EO gas sterilized

Code	Size
MEREN1607	0.16 x 7 mm
MEREN1613	0.16 x 13 mm
MEREN2007	0.20 x 7 mm
MEREN2013	0.20 x 13 mm
MEREN2207	0.22 x 7 mm
MEREN2213	0.22 x 13 mm

DongBang Acupuncture Needles

DongBang DB100



Silicone Coated Needles with Tubes & Spring Handles

- One needle with a guide tube per blister
- The needle features a spring handle and silicone coating
- Packed in 100 pcs per box
- EO gas sterilized - 5 year shelf life

Code	Size
DB1001615	0.16 x 15 mm
DB1001630	0.16 x 30 mm
DB1001815	0.18 x 15 mm
DB1001830	0.18 x 30 mm
DB1002015	0.20 x 15 mm
DB1002025	0.20 x 25 mm
DB1002030	0.20 x 30 mm
DB1002040	0.20 x 40 mm
DB1002050	0.20 x 50 mm
DB1002230	0.22 x 30 mm
DB1002240	0.22 x 40 mm
DB1002525	0.25 x 25 mm
DB1002530	0.25 x 30 mm
DB1002540	0.25 x 40 mm
DB1002550	0.25 x 50 mm
DB1003030	0.30 x 30 mm
DB1003040	0.30 x 40 mm
DB1003050	0.30 x 50 mm
DB1003540	0.35 x 40 mm



Acupuncture Needles

DongBang DB102



Silicone Coated Needles without Guide Tubes

- One needle without a guide tube per blister
- The needle features a spring handle and silicone coating
- Packed in 100 pcs per box
- EO gas sterilized - 5 year shelf life

Code	Size
DB1021615	0.16 x 15 mm
DB1021630	0.16 x 30 mm
DB1021815	0.18 x 15 mm
DB1021830	0.18 x 30 mm
DB1022015	0.20 x 15 mm
DB1022025	0.20 x 25 mm
DB1022030	0.20 x 30 mm
DB1022040	0.20 x 40 mm
DB1022050	0.20 x 50 mm
DB1022230	0.22 x 30 mm
DB1022240	0.22 x 40 mm
DB1022515	0.25 x 15 mm
DB1022525	0.25 x 25 mm
DB1022530	0.25 x 30 mm
DB1022540	0.25 x 40 mm
DB1023015	0.30 x 15 mm
DB1023030	0.30 x 30 mm
DB1023040	0.30 x 40 mm
DB1023050	0.30 x 50 mm
DB1023060	0.30 x 60 mm



DongBang Acupuncture Needles

DongBang DB105



Made from surgical stainless steel and EO gas sterilised, these sterile needles have a spring handle and easy-glide insertion action. The needles are gently released from the guide tube by removing the plastic tab.

- 5 needles and one guide tube per blister
- Disposable for single use
- 500pcs per box
- Silicon coated
- EO gas sterilized - 5 year shelf life

Code	Size
DB1051615	0.16 x 15 mm
DB1051630	0.16 x 30 mm
DB1051815	0.18 x 15 mm
DB1051830	0.18 x 30 mm
DB1052015	0.20 x 15 mm
DB1052025	0.20 x 25 mm
DB1052030	0.20 x 30 mm
DB1052040	0.20 x 40 mm
DB1052050	0.20 x 50 mm
DB1052230	0.22 x 30 mm
DB1052240	0.22 x 40 mm
DB1052515	0.25 x 15 mm
DB1052525	0.25 x 25 mm
DB1052530	0.25 x 30 mm
DB1052540	0.25 x 40 mm
DB1052550	0.25 x 50 mm
DB1053030	0.30 x 30 mm
DB1053040	0.30 x 40 mm
DB1053050	0.30 x 50 mm
DB1053060	0.30 x 60 mm



Acupuncture Needles

DongBang DB106



Made from surgical stainless steel and EO gas sterilised, these sterile needles have a spring handle and easy-glide insertion action. The needles are gently released from the guide tube by removing the plastic tab.

- 10 needles and one guide tube per blister
- Disposable for single use
- 1000pcs per box
- Silicon coated
- EO gas sterilized - 5 year shelf life

Code	Size
DB1061615	0.16 x 15 mm
DB1061630	0.16 x 30 mm
DB1061815	0.18 x 15 mm
DB1061830	0.18 x 30 mm
DB1061840	0.18 x 40 mm
DB1062015	0.20 x 15 mm
DB1062025	0.20 x 25 mm
DB1062030	0.20 x 30 mm
DB1062040	0.20 x 40 mm
DB1062050	0.20 x 50 mm
DB1062230	0.22 x 30 mm
DB1062240	0.22 x 40 mm
DB1062515	0.25 x 15 mm
DB1062525	0.25 x 25 mm
DB1062530	0.25 x 30 mm
DB1062540	0.25 x 40 mm
DB1062550	0.25 x 50 mm
DB1062560	0.25 x 60 mm
DB1063015	0.30 x 15 mm
DB1063030	0.30 x 30 mm
DB1063040	0.30 x 40 mm
DB1063050	0.30 x 50 mm
DB1063060	0.30 x 60 mm
DB1063530	0.35 x 30 mm
DB1063540	0.35 x 40 mm



DongBang Long Needles

- Long needle made from surgical stainless steel
- Ten spring handled needles with one tube per pouch



Code	Size	Contents
DB111105	0.35 x 105 mm	1000 pcs / box
DB111120	0.35 x 120 mm	1000 pcs / box
DB111135	0.35 x 135 mm	1000 pcs / box

DongBang Press Needles

- Press needles for ear and body acupuncture
- Made from surgical stainless steel



Code	Size	Contents
DB130S	0.20 x 1.0 mm	100 pcs / box
DB130M	0.20 x 1.5 mm	100 pcs / box

DongBang Intradermal Needles

- Silicone-free, stainless steel intradermal needles
- Packed in 50 pcs per box and 5 pcs per pouch



Code	Size
DB131S	0.12 x 3 mm
DB131M	0.12 x 5 mm
DB131L	0.14 x 7 mm

DongBang Hand Needles

- Hand acupuncture needles featuring spring handles
- Packed in 100 pcs per box and 10 pcs per pouch



Code	Size
DB132	0.18 x 8 mm

Acupuncture Needles

ASP Original Needles

ASP Original Auriculotherapy Needles



The ASP Original needle is composed of a small stainless steel needle housed in a plastic injector. This injector allows the practitioner to easily and precisely implant the needle into the patient's ear. Inserting the needle in this way avoids any risk of contamination. After placement, the needle can be protected by an adhesive plaster. The needle ensures a permanent action for as long as it remains in place. A black magnet housed in the injector can be turned around in front of the implanted needle. ASP Original needles are gamma ray sterilized and for single use only.

Invented in collaboration with Dr Paul Nogier, the ASP Original or "Aiguille Semi-Permanente" is used to perform ear acupuncture. It has been marketed for more than 30 years by Sedatelec, and the shape of this small needle makes it exceptionally efficient.

Code	Size
ASC8	Classic 8 pcs/box
ASC80	Classic 80 pcs/box
ASC200	Classic 200 pcs/box
ASG8	Gold 8 pcs/box
ASG80	Gold 80 pcs/box
ASG200	Gold 200 pcs/box

ORCA™ Acupuncture Needles

ORCA™ Accent Copper Handled Needles

- Chinese-style needle with copper wire handle
- One needle with a tube per blister
- Packed in 100 pcs per box



Code	Size
OCT2525	0.25 x 25 mm
OCT2530	0.25 x 30 mm
OCT3030	0.30 x 30 mm
OCT3040	0.30 x 40 mm
OCT3050	0.30 x 50 mm
OCT3075	0.30 x 75 mm

ORCA™ Advanced Copper Handled Needles

- Chinese-style needle with copper wire handle
- One needle without tube per blister
- Packed in 100 pcs per box



Code	Size
OC2015	0.20 x 15 mm
OC2515	0.25 x 15 mm
OC2525	0.25 x 25 mm
OC2540	0.25 x 40 mm
OC3030	0.30 x 30 mm

ORCA™ Accent Silver Handled Needles

- Chinese-style needle with silver-plated wire handle
- One needle with a tube per blister
- Packed in 100 pcs per box



Code	Size
OST1615	0.16 x 15 mm
OST1630	0.16 x 30 mm
OST2015	0.20 x 15 mm
OST2525	0.25 x 25 mm
OST2530	0.25 x 30 mm
OST2540	0.25 x 40 mm
OST3030	0.30 x 30 mm
OST3040	0.30 x 40 mm
OST3050	0.30 x 50 mm
OST3060	0.30 x 60 mm
OST3070	0.30 x 70 mm
OST3075	0.30 x 75 mm
OST3080	0.30 x 80 mm
OST3090	0.30 x 90 mm
OST30100	0.30 x 100 mm

ORCA™ Advanced Silver Handled Needles



Chinese-style needles with silver-plated wire handles

- One needle without guide tube per blister
- Packed in 100 pcs per box
- Features silicone coating for easy insertion
- EO gas sterilized

Code	Size
OS2015	0.20 x 15 mm
OS2525	0.25 x 25 mm
OS2530	0.25 x 30 mm
OS2540	0.25 x 40 mm
OS3030	0.30 x 30 mm
OS3040	0.30 x 40 mm

ORCA™ Natural Copper Handled Needles



Silicone-free needles with copper wire handles

- One needle without guide tube per blister
- Packed in 100 pcs per box
- Ideal for electro-acupuncture
- EO gas sterilized

Code	Size
ON2525	0.25 x 25 mm

ORCA™ Acupuncture Needles

ORCA™ Resolve Silver Handled Needles



Chinese-style needles with silver-plated wire handles

- Five needles with one guide tube per blister
- Packed in 500 pcs per box
- Features silicone coating for easy insertion
- EO gas sterilized

Code	Size
OS52525	0.25 x 25 mm
OS52530	0.25 x 30 mm

ORCA™ Resolve Copper Handled Needles



Chinese-style needles with copper wire handles

- Five needles with one guide tube per blister
- Packed in 500 pcs per box
- Features silicone coating for easy insertion
- EO gas sterilized

Code	Size
OC52525	0.25 x 25 mm
OC52530	0.25 x 30 mm

ORCA™ Premium Plastic Handled Needles



Coloured-coded plastic handled needles

- One needle without tube per blister
- Packed in 100 pcs per box

Code	Size
OP2015	0.20 x 15 mm
OP2525	0.25 x 25 mm
OP2540	0.25 x 40 mm
OP3030	0.30 x 30 mm

Acupuncture Needles

Seirin® Acupuncture Needles

Seirin® J-Type Needles



- One needle with a guide tube per blister
- Packed in 100 pcs per box
- Every needle features a skin-friendly plastic handle
- EO gas sterilized

Code	Size
SJ1215	0.12 x 15 mm
SJ1230	0.12 x 30 mm
SJ1415	0.14 x 15 mm
SJ1430	0.14 x 30 mm
SJ1615	0.16 x 15 mm
SJ1630	0.16 x 30 mm
SJ1830	0.18 x 30 mm
SJ1840	0.18 x 40 mm
SJ2030	0.20 x 30 mm
SJ2040	0.20 x 40 mm
SJ2050	0.20 x 50 mm
SJ2530	0.25 x 30 mm
SJ2540	0.25 x 40 mm
SJ2550	0.25 x 50 mm
SJ3030	0.30 x 30 mm
SJ3040	0.30 x 40 mm
SJ3050	0.30 x 50 mm



Seirin® Spinex Intradermal Needles



- One EO gas sterilized needle per blister
- Packed in 100 pcs per box
- Features a ring handle measuring 1.5 mm in diameter
- Made from surgical stainless steel

Code	Size
DBSSP1203	0.12 x 3 mm
DBSSP1204	0.12 x 4 mm
DBSSP1406	0.14 x 6 mm

Seirin® B-Type Needles



- One EO gas sterilized needle per blister
- Packed in 100 pcs per box
- Features silicone coating for easy insertion
- Made from surgical stainless steel
- Every B-Type needle features a lightweight plastic handle

Code	Size
SB1615	0.16 x 15 mm
SB2015	0.20 x 15 mm
SB2540	0.25 x 40 mm
SB3030	0.30 x 30 mm

Seirin® New Pyonex Needles



The patented plastic packaging allow the practitioner to insert the needle without touching the needle tip or the surface of the Micropore™ adhesive tape.

The skin friendly adhesive allows the skin to breathe during treatment. With Seirin's smooth insertion, the result is a painless experience. Packed in 100 pcs / box.

Code	Size	Colour
DBSNP2003	0.20 x 0.3 mm	Orange
DBSNP2006	0.20 x 0.6 mm	Yellow
DBSNP2009	0.20 x 0.9 mm	Green
DBSNP2012	0.20 x 1.2 mm	Blue
DBSNP2015	0.20 x 1.5 mm	Pink

Acupuncture Needles

Acupuncture Accessories

Hand Needle Injectors



- Assists in the insertion of hand needles
- Magnets inside pick up, hold and secure the needle
- The weight of the plunger will drop, using gravity or the spring force to insert the needle
- DB135 allows adjustments for the depth of insertion

Code	Type
DB134	Gravity
DB135	Spring Force

NIDO 7-Star Needles



- Disposable 7-star needle
- Pipe type double-ended – ideal for bleeding acupuncture
- Individually packed 7-star dermal needle
- Double-side needle-18cm(7") long with 7 needles

Code	Type
ND151	1 head

Lancet Devices

DongBang NANOLET™ Lancets



Code	Size
DB40023G	0.60 mm (23G)
DB40026G	0.45 mm (26G)
DB40028G	0.35 mm (28G)
DB40030G	0.30 mm (30G)

LANZO



Adjustable tip - 5 settings for custom lancet penetration
Lancet ejector - no need to touch the lancet to remove it
Reduced pain - thanks to the linear tracking system
Silent - the carrier is provided with two springs and a stopper

Optional clear tip for alternative site testing designed in order to reduce pain and secure a blood sample.

DB403A

DongBang Lancet Device



- To be used with DB400 NANOLET™ Lancets
- Press the trigger on the top of the pen to push the needle

Code
DB403B

Acupuncture Needles

Park Sham Device

Park Sham Research Tools



Fig 1


The Park Sham Device is an apparatus used to administer the sham needle and real needles. The Park Sham Needle appears similar to the real needle, but is blunt and free to slide within the handle so that when pressed, it telescopes into the handle instead of penetrating the skin. The device holds the un-inserted sham needle in place and maintains the sterility of a genuine acupuncture needle.

This device consists of a guide tube: a special oversized tube with base flange which is equipped with double-sided tape and can therefore stick on target skin (Fig 1). Since the Park Sham Device was first developed in 1999, several teams of researchers in UK, Japan, Korea and USA have tested its credibility and validity in various settings.

Based on these findings over the past 19 years it has been frequently used in research projects, often funded by government agencies including the US National Institute of Health. Park and his fellow colleagues, who developed this device, state that this device is indistinguishable from the same procedure using real needles in acupuncture naive subjects and inactive, where the specific needle sensation is taken as a surrogate measure of activity. However, due to the nature of minimal physiological activity required to achieve successful subject blinding in acupuncture trial, they emphatically call for caution in application of sham acupuncture control and in interpretation of the findings. For this reason, Park is keen to call this a bona fide version of sham control rather than final version.

- Park J. et al. Development of a new sham needle, *Acupunct Med* 17 (1999), pp. 110-112.
- Park J. et al., Validating a new nonpenetrating sham acupuncture device: two randomised controlled trials, *Acupunct Med* 20 (2002), pp. 168-174.
- Whale CA et al. Pilot study to assess the credibility of acupuncture in acute exacerbations of chronic obstructive pulmonary disease. *Acupunct Med*. 2009 Mar;27(1):13-5.
- Tsukayama H et al. Factors that influence the applicability of sham needle in acupuncture trials: two randomized, single-blind, crossover trials with acupuncture-experienced subjects. *Clin J Pain*. 2006 May;22(4):346-9.
- Chae Yet al. The neural substrates of verum acupuncture compared to non-penetrating placebo needle: an fMRI study. *Neurosci Lett*. 2009 Jan 30;450(2):80-4.
- Tan CW et al. Discrimination of real and sham acupuncture needles using the Park sham device: a preliminary study, *Arch Phys Med Rehabil* 90 (2009), pp. 2141-2145.
- Downs NM et al. The effect of real and sham acupuncture on thermal sensation and thermal pain thresholds. *Arch Phys Med Rehabil* 2005;86:1252-7.
- Park J, et al. Acupuncture for subacute stroke rehabilitation: a Sham-controlled, subject and assessor blind, randomized trial. *Arch Intern Med*. 2005 Sep 26;165(17):2026-31.
- Smith Pet al. The efficacy of acupuncture in the treatment of temporomandibular joint myofascial pain: a randomised controlled trial. *J Dent*. 2007 Mar;35(3):259-67.
- Kennedy Set al. Acupuncture for acute non-specific low back pain: a pilot randomised non-penetrating sham controlled trial. *Complement Ther Med*. 2008 Jun;16(3):139-46.
- Park J. Sham needle control needs careful approach. *Pain*. 2004 May;109(1-2): 195-6
- Park JJ. Developing and validating a sham acupuncture needle. *Acupunct Med*. 2009 Sep;27 (3):93.

How to use the Park Sham Device



Code	Product
PSD01	Park Sham Device
PSD03	Sham Needle Only
PSD04	Sticky Tape
PSD05	Tube Only
PSD06	Real Matching Needle

To avoid compromising the sterility, you could peel off only the top of the package to expose the handle of the needle.

Peel off the cover of the double-sided sticky tape and place the device on the skin while keeping the needle handle projecting about 5mm from the top of the Park tube.

Then tap the needle gently to insert it about 5mm followed by moving down the upper telescopic guide tube to allow sufficient exposure of the handle for manipulation of the needle.

To remove: first take out the needle and dispose of it safely, then peel the flange from the skin, with the guide tube in place. Doing this quickly will minimise discomfort.


Kinesiology Tape Working Principles



Kinesiology taping was developed in the 1970s by Japanese chiropractor Kenzo Kase, who went to the U.S. in the mid-1990s. Unlike standard athletic taping, which often involves wrapping a joint for support and compression, kinesiology tape is placed in a variety of patterns depending on the injury. It is pulled to differing degrees of tension to create the desired effect and is typically worn for two to five days, unlike standard tape, which is used mainly during an activity.

Kinesiology tape can be used to pull back a shoulder that is hunching forward or to reduce swelling in a joint. Tape can be used to pull up the skin and create an area of low pressure where fluid can move and drain. Some clinicians believe the tape provides stimulation to skin cells that affects pain pathways - similar to rubbing a spot that hurts. (Summarized from the Wall Street Journal Apr 27, 2010)

Directions for use



Always shape the cut ends of the tape in such a way that they are curved with no sharp corners. This helps to reduce the problem of the corners catching in clothing, causing the tape to become dislodged.

To correctly apply the tape, first measure the required length while the backing tape is still in place, then shape the corners as outlined above. Position the area to be taped so that the underlying muscles are fully stretched. Remove the backing and lay the tape on the skin sticking it in place. Do not stretch the tape, simply lay it on the skin.

The exact length of time it remains in place will depend on how greasy the skin is and how much the subject perspires. It is possible to shower carefully, keeping the tape in place. The area should be patted dry and not rubbed. If the tape then feels unpleasantly damp, the area may be dried using a hairdryer on the lowest heat setting, held some distance from the body.

Kinesiology Tapes

meridius® ADAPTA™ Tape



Meridius® Adapta Tape is the result of years of development and testing conducted by top athletes and renowned scientists. This kinesiology tape is comfortable to wear, but also sticks for up to 4 days due to the premium hypo-allergenic adhesive that's being used.

Using only high quality materials, we have created a water resistant, bi-linearly elastic tape that doesn't compromise the natural range of motion. Adapta Tape has shown results in improving circulation, helping with muscle soreness, supporting joints and improving posture.

- Desert Beige
- Chaos Black
- Arctic Blue
- Toxic Pink
- Mystic Violet
- Atomic Yellow
- Viper Green
- Tiger Orange
- Volcano Red

Code	Size	Colour
ME01BE	5 m	Desert Beige
ME03BK	5 m	Chaos Black
ME05BL	5 m	Arctic Blue
ME07PK	5 m	Toxic Pink
ME09YL	5 m	Atomic Yellow
ME11GR	5 m	Viper Green
ME12OR	5 m	Tiger Orange
ME13RD	5 m	Volcano Red
ME14VT	5 m	Mystic Violet

Code	Size	Colour
ME32BE	32 m	Desert Beige
ME32BK	32 m	Chaos Black
ME32BL	32 m	Arctic Blue
ME32PK	32 m	Toxic Pink
ME32YL	32 m	Atomic Yellow

AcuTape™



This tape is made of:

- Cotton 72%
- Spandex 3%
- Ethyl acetate 15%

Code	Size	Colour
DAT01BE	5 m	Beige
DAT03BK	5 m	Black
DAT05BL	5 m	Blue
DAT07PK	5 m	Pink
DAT09YL	5 m	Yellow
DAT11GR	5 m	Green

meridius® KROS™ Tape



This lattice-shaped, adhesive-backed therapeutic patch works by targeting myofascial trigger-points, or "knots" in the muscular tissues, which often cause far-reaching pain and discomfort.

By applying constant, gentle pressure to specific points, KROS™ tape helps encourage blood flow, reduce pain and enhance proprioception.

- Therapeutic pressure to trigger points
- Reduce pain and encourage blood flow
- Long lasting adhesive stays put for days
- Cross hatch design allows breathability
- Water resistant
- Wraps around curves and joints

Code	Size
ME301S	Small, 180 pcs
ME301M	Medium, 120 pcs
ME301L	Large, 40 pcs


Non-Stick Scissors



- Teflon-coated scissors
- Ideal for cutting and shaping kinesiology tapes
- Developed to minimise glue build-up on the blades

Code	Type
AT30	Teflon-coated, 5"

Acupressure Taping Book



Learn about the combination of acupuncture and taping in this comprehensive book on the practice of acupressure taping. Authors Hecker and Liebchen offer in-depth advice and practical guidance on acupressure taping through full colour step-by-step guides and detailed descriptions. From preventative taping to therapeutic techniques, Acupressure Taping is one of the best books available on the subject.

AT20

meridius® Cohesive Bandage



- Ideal for compression, strapping and holding ice packs
- Comes in woven and non-woven versions

Code	Size	Colour / Type
MECA03-25BE	2.5 cm x 4.5 m	Beige / Standard
MECA03-50BE	5.0 cm x 4.5 m	Beige / Standard
MECA03-75BE	7.5 cm x 4.5 m	Beige / Standard
MECA02-25BE	2.5 cm x 4.5 m	Beige / Non Woven
MECA02-50BE	5.0 cm x 4.5 m	Beige / Non Woven
MECA02-50BK	5.0 cm x 4.5 m	Black / Non Woven
MECA02-50RD	5.0 cm x 4.5 m	Red / Non Woven
MECA02-75BE	7.5 cm x 4.5 m	Beige / Non Woven
MECA02-75BK	7.5 cm x 4.5 m	Black/ Non Woven
MECA02-75BL	7.5 cm x 4.5 m	Blue / Non Woven
MECA02-75RD	7.5 cm x 4.5 m	Red / Non Woven

meridius® Crepe Bandage



- Cohesive high compression bandage made of crepe
- Ideal to be used as an aid in the treatment of sprains
- Adapts to body contours making it comfortable to wear
- Made of highly elastic material

Code	Size
MECA04-10	10 cm x 4 m
MECA04-05	5 cm x 4 m

meridius® Elastic Adhesive Bandage



- Made of high quality elastic cotton fabric
- Ideal for ankle, knee, hand and wrist taping
- Features strong, hypoallergenic adhesive
- Can be used in rugby for a better grip on line-outs

Code	Size	Type
ME34EAB25	2.5 cm x 4.5 m	Standard EAB
ME34EAB50	5.0 cm x 4.5 m	Standard EAB
ME34EAB75	7.5 cm x 4.5 m	Standard EAB
METAB23	2.5 cm x 6.9 m	Tear EAB
METAB24	5.0 cm x 6.9 m	Tear EAB
METAB25	7.5 cm x 6.9 m	Tear EAB

meridius® Zinc Oxide Tape



- Ideal for prevention and treatment of ligament sprains
- Features strong, hypoallergenic adhesive
- Can be used to protect areas from blisters
- Commonly used by weightlifters

Code	Size
MECA05-25	2.5 cm x 10 m
MECA05-38	3.8 cm x 10 m
MECA05-50	5.0 cm x 10 m

SpoBand® Elastic Resistance Band



SpoBand® is a premium resistance band that is used for muscle toning, posture correction, cardio exercises, rehabilitation and yoga.

While the typical rubber or latex bands can begin to smell after usage, SpoBand® doesn't. It's also easy to clean as it is machine washable.

Another important feature is that SpoBand® doesn't stick and melt as the regular bands do. Moreover, the high quality materials will not tear or snap like latex bands are prone to.

SpoBand features 9 loops for different resistances making it versatile, safe and easy to use.

Code	Name	Size
MESB01	SpoBand 20	110 x 2.0 x 0.5 cm
MESB02	SpoBand 25	110 x 2.5 x 0.5 cm
MESB03	SpoBand 30	110 x 3.0 x 0.5 cm
MESB04	SpoBand 35	110 x 3.5 x 0.5 cm
MESB05	SpoBand 45	110 x 4.5 x 0.5 cm
MESB06	SpoBand 55	110 x 5.5 x 0.5 cm

meridius® Cupping Sets

meridius® 30pcs CuppingSet



Meridius® Cupping Sets are presented in three sizes: 10 pieces, 17 pieces and 30 pieces. They come in a blue or black case with a firm handle which is ideal for traveling and storage. The blue silicone valve can be detached from the cup for easy sterilisation and maintenance.

Our cupping sets are extremely versatile due to the large number of cupping shapes and dimensions available. Made from high-quality polycarbonate, Meridius® Cupping Sets are extremely durable and do not crack from falling or prolonged usage.

ME17 Contents

Code	Diameter	Pieces
ME17-RO1	Round, 55 mm	8
ME17-CO2	Contoured, 43 mm	1
ME17-CO3	Contoured, 41 mm	1
ME17-CO4	Contoured, 34 mm	1
ME17-RO4	Round, 27 mm	1
ME17-TR	Triangle, 23 mm	1
ME17-OV1	Contoured, 72 mm	1
ME17-CO1	Contoured, 61 mm	1
ME17-LO1	Round, 50 mm	1
ME17-LO2	Round, 40 mm	1
ME17-T	Extension Tube	1
ME17-P	Hand Pump	1

ME10 Contents

Code	Diameter	Pieces
ME17-RO1	Round, 55 mm	6
ME17-CO2	Contoured, 43 mm	2
ME17-CO1	Contoured, 61 mm	1
ME17-CO4	Contoured, 34 mm	1
ME17-P	Hand Pump	1

Cupping Therapy & Massage

meridius® Cupping Sets

meridius® 17pcs Cupping Set



meridius® 10pcs Cupping Set



Meridius® Cupping Sets are presented in three sizes: 10 pieces, 7 pieces, and 30 pieces. They come in a blue or black case with a firm handle which is ideal for traveling and storage. The blue silicone valve can be detached from the cup for easy sterilisation and maintenance.

Our cupping sets are extremely versatile due to the large number of cupping shapes and dimensions available. Made from high-quality polycarbonate, Meridius® Cupping Sets are extremely durable and do not crack from falling or prolonged usage.

ME17 Contents

Code	Diameter	Pieces
ME17-R01	Round, 55 mm	8
ME17-C02	Contoured, 43 mm	1
ME17-C03	Contoured, 41 mm	1
ME17-C04	Contoured, 34 mm	1
ME17-R04	Round, 27 mm	1
ME17-TR	Triangle, 23 mm	1
ME17-OV1	Contoured, 72 mm	1
ME17-C01	Contoured, 61 mm	1
ME17-L01	Round, 50 mm	1
ME17-L02	Round, 40 mm	1
ME17-T	Extension Tube	1
ME17-P	Hand Pump	1

ME10 Contents

Code	Diameter	Pieces
ME17-R01	Round, 55 mm	6
ME17-C02	Contoured, 43 mm	2
ME17-C01	Contoured, 61 mm	1
ME17-C04	Contoured, 34 mm	1
ME17-P	Hand Pump	1

Cupping Therapy & Massage

Replacement cups



ME17-RO1



ME17-CO2



ME17-CO3



ME17-CO4



ME17-RO4



ME17-TR



ME17-OV1



ME17-CO1



ME17-LO1



ME17-LO2



ME17T



ME17P

Code	Shape	Diameter
ME17-RO1	Round	Round, 55 mm
ME17-CO2	Contoured	Contoured, 43 mm
ME17-CO3	Contoured	Contoured, 41 mm
ME17-CO4	Contoured	Contoured, 34 mm
ME17-RO4	Round	Round, 27 mm
ME17-TR	Triangle	Triangle, 23 mm
ME17-OV1	Contoured	Contoured, 72 mm
ME17-CO1	Contoured	Contoured, 61 mm
ME17-LO1	Round	Round, 50 mm
ME17-LO2	Round	Round, 40 mm
ME17-T	Extension Tube	Standard
ME17-P	Hand Pump	Standard

DongBang Cupping Sets

DongBang Plastic Cupping Set



- Elegant lightweight case - Easy to carry and store
- Cups are non-breakable with valve on the top
- Place the pump on top of the cup and apply pressure
- To release pressure press the valve

Code	Diameter
DB350	51 mm x 13 pcs
	46 mm x 1 pcs
	41 mm x 1 pcs
	33 mm x 1 pcs
	25 mm x 1 pcs
	Tube + Hand Pump
DB360	51 mm x 6 pcs
	46 mm x 1 pcs
	41 mm x 1 pcs
	33 mm x 1 pcs
	25 mm x 1 pcs
	Hand Pump

DongBang Disposable Cups



- Disposable plastic cups for single use
- Compatible with DB304 hand pistol pump

Code	Diameter	Contents
DB310-1	55 mm	5 cups / pack
DB310-2	51 mm	5 cups / pack
DB310-3	46 mm	5 cups / pack
DB310-4	41 mm	5 cups / pack
DB310-5	36 mm	5 cups / pack

Cupping Therapy & Massage

Bio Magnet Cupping Set



- A 17 piece plastic cupping set with magnets which also apply a mild acu-pressure to the body
- Each set includes 7 red and 8 blue colour coded, 0.3" dia, 1,000 gauss, spring adjustable, 0.3" dia x 1.7" long, removable magnetic insert probes. These fit the 2" and 1.5" dia cups
- Also includes pistol grip hand pump, extension tube and carrying case

Code
DBA341

Plastic Cupping Set



Code	Size
DBA340-6	65 mm x 1 pcs
	55 mm x 1 pcs
	45 mm x 1 pcs
	35 mm x 1 pcs
	Tube + Pump
DBA340-12	65 mm x 2 pcs
	55 mm x 2 pcs
	45 mm x 2 pcs
	35 mm x 4 pcs
	Tube + Pump
DBA340-24	65 mm x 12 pcs
	55 mm x 4 pcs
	45 mm x 4 pcs
	35 mm x 4 pcs
	Oil + Tube + Pump

Glass Cupping Sets

DaeGun Glass Cupping Set



- Silver colour carrying case with prime quality glass cups
- Offers better cupping treatment
- Glass cup with valve on the top
- Contains 14 cups with a hand pump

DB302

Glass Fire Cupping Jars



- Thick, heavy glass jar for cupping therapy
- Ideal for fire cupping
- Burn alcohol cotton balls to create a vacuum
- Easy to clean

Code	Diameter
DB305-XL	75 mm
DB305-L	65 mm
DB305-M	55 mm
DB305-S	45 mm
DB305-XS	40 mm

Cupping Therapy & Massage

Other Cupping Sets

HACI® Magnetic Cupping Set



- Rubber bulb plastic suction cups
- 2300 gauss strength magnets
- 5 bio-north cups and 5 bio-south cups

Code
DBA320

Silicone Cupping Jars



- Each set contains 4 cups (25, 40, 50, 60mm diameter)

Code	Colour
DBA322T	Transparent
DBA322	Blue

Bamboo Cupping Jars



- Natural bamboo cups in three different sizes
- Cups can be soaked in water for infusion with herbs

Code	Contents
DBA321	3 pcs / box

Silicone Cups for Face and Eyes



Code	Size	Colour
DBA325S-BL	6 mm x 24 mm	Blue
DBA325S-GR	6 mm x 24 mm	Green
DBA325S-PK	6 mm x 24 mm	Pink
DBA325S-PL	6 mm x 24 mm	Purple
DBA325S-CL	6 mm x 24 mm	Clear

Cupping Therapy & Massage

Silicone Cups for General Use



Code	Size	Colour
DBA324L-01BL	47 mm x 80 mm	Blue
DBA324L-01CL	47 mm x 80 mm	Clear
DBA324L-01GR	47 mm x 80 mm	Green
DBA324L-01PK	47 mm x 80 mm	Pink
DBA324L-01PL	47 mm x 80 mm	Purple
DBA324M-01BL	33 mm x 60 mm	Blue
DBA324M-01CL	33 mm x 60 mm	Clear
DBA324M-01GR	33 mm x 60 mm	Green
DBA324M-01PK	33 mm x 60 mm	Pink
DBA324M-01PL	33 mm x 60 mm	Purple
DBA324S-01BL	19 mm x 80 mm	Blue
DBA324S-01CL	19 mm x 80 mm	Clear
DBA324S-01GR	19 mm x 80 mm	Green
DBA324S-01PK	19 mm x 80 mm	Pink
DBA324S-01PL	19 mm x 80 mm	Purple

Medical grade silicone, odour free and hygienic. Transparency enables observation of the skin and soft tissue inside the cups during the application.

How do they work?

Expelling air out of the cup creates negative pressure and causes the skin and soft tissue to rise up. Thanks to the soft rims of the silicone cups, these cups can be applied to contoured and other tricky areas of the body. The soft rim creates a flawless seal around the connected skin and is far more ergonomic than its glass counterparts.

How are they used?

Ensure the skin is clean. Apply massage cream, gel or oil to moisten the skin and prepare it for both massaging and cup gliding. Squeeze the cup in the middle to release air from inside. Place the squeezed cup onto the skin and hold it against the skin to ensure a complete seal is created.

Slowly release and allow the pressure to lift the skin and flesh. Grasp the base of the cup and gently massage or glide along the skin to create a wave of lifting. Massage oil, cream, gel, or soapy water can be considered to keep the skin lubricated to assist the massaging.

Non-Sterile Cotton Wool Balls



Code	Size
DBA6059000	Small, 500 pcs / pack
DBA6059001	Large, 250 pcs / pack

Forceps



Code	Size
DBA60518033	25 cm / 10"
DBA605763	18 cm / 7"

Iris 10cm Non-Toothed Forceps



Code	Size
DBA605716	10 cm / 4"

12cm Adson Non Toothed Forceps



Code	Size
DBA605718	12 cm / 4 ³ / ₄ "

Gua Sha Board



- Flat rectangle with slight contour at sides for traditional Gua Sha treatments

Code	Size
BN6007	3.5 x 6.5 cm
BN6008	4.5 x 8 cm
BN6009	5.5 x 10 cm

Gua Sha Jade



Code	Size
BN6011H	5 x 7.5 cm
BN6011C	5.5 x 8 cm
BN6011F	4 x 11 cm

Gua Sha Rod



Code
BN6012

Gua Sha Roller Jade



Code
BN6015D

Massage Oil with Soya Bean & Wheatgerm Oil



- Light yellow colour. Available in 500ml or 4L
- Contains: soya bean and wheatgerm oils

Code	Size
HT6027	500 ml
HT6027	4000 ml

Press Pellet & Probes

Press Pellet

Magnetic Ball Plasters



- 800 gauss strength
- 2mm diameter with tan-coloured surgical tape backing

Code	Type	Contents
DBA424	Silver plated, 2mm	100pcs / box
DBA425	Gold plated, 2mm	100pcs / box

Magrain Ion Pellet



- Sakamura Magrain are used for both ear and body points
- Magrain are silver colour 1.2mm alloy granules
- Magrain is unique in that it provides an ionizing reaction which further enhances the compression stimulation to the point

Code	Type	Contents
DB509006C	Non-plated	300pcs / pack
DB509006GC	Gold plated	300pcs / pack

Vaccaria Ear Seeds



- Vaccaria Ear Seeds are attached to an individual plaster, so that after the removal they can be placed immediately
- Proper application on the desired place is possible thanks to the plaster
- Typical application of the ear seed is on the acupuncture points of the ear, in order to achieve the lasting effect of the pressure
- Based on their size they can also be attached to other locations of the body

Code	Size	Contents
DBA426	2 mm	100pcs / box

Press Pellet



- Ideal for Korean Hand Acupuncture Therapy
- The metal protrusions provide gentle stimulation
- Pellets may remain in place for several hours

Code	Type	Contents
DB506A	Silver plated	100pcs / box
DB506B	Silver plated	30pcs / box
DB506C	Gold plated	100pcs / box
DB506D	Gold plated	30pcs / box

Tourmaline Stone Granules for Ear



- The product must be used by a qualified practitioner
- Product is for single use only, with care, to attach securely to ear
- Store in a safe place away from children, chemicals, moisture and extreme temperature

Code	Type	Contents
DBA531	1 ball	420pcs / box
DBA531-123	1 ball	126pcs / box
	2 balls	84pcs / box
	3 balls	140pcs / box

Press Pellet & Probes

Magnets


HANIL Magnetic Press Pellet



- Magnets for acupuncture application
- 5mm diameter pellet on plaster backing
- 800 gauss
- Contains 100 pcs per box

Code	Diameter	Contents
DB511	5 mm	100pcs / box

ACCU-BAND 800G



ACCU-BAND gives a constant pressure and a concentrated magnetic flux of 800 gauss, 6000 gauss, 9000 gauss. Adheres even in water and will not rust. The south pole faces towards the body.

Code	Type
DB509001G	Gold plated, 800 gauss, 24pcs/ box
DB509001S	Non-plated, 800 gauss, 24pcs/ box
DB509002S	Gold plated, 6000 gauss, 12pcs/ box
DB509003S	Gold plated, 9000 gauss, 12pcs/ box

Hpro Ferrite Magnets



Code	Type
DBA5901	Non-plate, 700 gauss 10pcs/Box
DBA5901G	Gold-plated, 700 gauss 10pcs/Box
DBA5902	Non-plate, 1000 gauss 10pcs/Box
DBA5902G	Gold-plated, 1000 gauss 10pcs/Box
DBA5905	Non-plated, 1200 gauss 10pcs/Box
DBA5905G	Gold-plated, 1200 gauss 10pcs/Box
DB5903G	Gold-plated, 2500 gauss 10pcs/Box

Replacement Plasters



- Small (13 mm) is used for all kinds of press pellet
- Medium (19 mm) is used for the magnetic press needles
- Large (30mm) was created for use with the 6,000 gauss magnetic (MEA) electrodes or rubber GSR electrodes

Code	Diameter	Contents
IT5007S	13 mm	96pcs / pack
IT5007M	19 mm	120pcs / pack
IT5007L	30 mm	80pcs / pack

Magnetic Needle 3000 Gauss



- Magnetic needle, able to reach most of the auricular and body points
- The other end is rounded and may also be used for acu-pressure and acu-massage

BN52500	3000 gauss	9 x 1.1 cm
---------	------------	------------

Magnetic Massage Balls



- Chrome-plated plastic balls
- Spiky massage balls for hands and body

Code	Type
HT6019	1000 gauss

Spring Loaded Dermal Roller with Probe



- Multi-functional, dermal roller with spring active probe
- Rolling heads with pyramid-shaped surface

Code	Size
DB502B	Handle: 50 x 13 mm Roller: 25 x 13 mm

Regular Dermal Roller



- Stainless steel
- Longer round handle

Code	Size
DB507	Rolling head: 13 x 25 mm

Su Jok Probes

Su Jok Diagnostic Probe



- Light and convenient probe, made of polystyrene

Code	Size
DBA5300	125 mm

Sharp Su Jok Probe



- Made of Chromium-nickel steel

Code	Size
DBA5302	100 x 1.5 x 2.5 mm

Round Su Jok Probe




- Made of Chromium-nickel steel

Code	Size
DBA5301	100 x 2.5 x 4.5 mm

Press Pellet & Probes

Su Jok Probes

Pressure Feelers



- To detect painful points of the ear
- Four different weights to ascertain the sensitivity of the irritated point.

Code	Colour
AS500BK	Black, 130g
AS500BL	Blue, 250g
AS500OR	Orange, 400g
AS500WT	White, 650g

Detector



- To locate and determine energy value of points
- Points detected with black, gold or north need to be treated with gold
- Points detected with white, silver or south need to be treated with silver

Code	Colour
AS501BW	Black / White
AS501GS	Gold / Silver
AS501NS	North / South

Ear Glass Rod



- Stimulates ear points by using glass rods
- 2 pcs per set

Code
AS502

Long Single Probe

- Very long, single probe for easy grip and application
- The handle is chiselled and offers good grip



Code

DB501A

Dual Probe

- The Dual-Probe includes pyramid spikes in the centre of the shaft.



Code

DB501B

Spring Loaded Probe

- Popular spring activated cutaneous probe
- Ideal for treating or locating trigger points



Code

DB502A

Tsun Meter Point Maker

- Interprets 1 Tsun into different distances to easily locate acupuncture points



Code

MM5001


Finger Massage Rings



Stainless steel and elastic enable the massage of fingers. Since a projection of the entire body is on each finger, the massage with a massage ring has a tangible strengthening and healing effect.

Code	Size	Colour
DBA5307BL	27 x 0.5 mm	Blue
DBA5307GR	27 x 0.5 mm	Green
DBA5307OR	27 x 0.5 mm	Orange
DBA5307RD	27 x 0.5 mm	Red
DBA5307YW	27 x 0.5 mm	Yellow
DBA5307WT	27 x 0.5 mm	White
DBA5307GD	27 x 0.5 mm	Gold

Massage Balls with Finger Ring



The massage balls are made of polystyrene. Massages hands by rolling/squeezing in a fist. Each set comes with 2 massage rings.

Code	Size	Colour
DBA5306BL	35 mm	Blue
DBA5306GR	35 mm	Green
DBA5306RD	35 mm	Red
DBA5306YL	35 mm	Yellow

Hand Massage Roller



Massage bracelet made from a resilient stainless steel wire. Used for massaging the whole hand and the wrist. The bracelets are elastic and suitable for hands of any sizes.

Code	Colour
DBA5308BL	Blue
DBA5308GR	Green
DBA5308OR	Orange
DBA5308RD	Red
DBA5308YW	Yellow
DBA5308WT	White
DBA5308GD	Gold
DBA5308SL	Silver

Small Spiral Massager



Stainless steel wire, used for the massage of both sides of the hands. The spirals are flexible and follow the hand lines, reaching even the most challenging areas.

Code	Size
DBA5305S	100 x 14 mm
DBA5305M	80 x 8.5 mm

Palm-Sized Three Channel Electric Stimulator



- For busy clinics that require several units
- Palm-sized and lightweight
- Three output channels
- Low/High intensity settings
- Unique asymmetric bi-phasic pulse
- Includes 3 sets of crocodile, needle connecting wires

Features

- High performance and versatility
- Low/High intensity settings
- Frequency selector offers a complete frequency range 1 to 500 Hz
- Suitable for needles or magnets
- Unique wave form
- Unique wave pulse minimises excessive charge build up in the tissue and damaging the nerve
- Each of the 3 channels are independently controlled
- Special dial controls allow easy use yet safe operation
- Built in battery life check with indicator light

Code
IT7002

Six-Channel Electroacupuncture Unit



- Power source: Four alkaline 'LR14' cells
- Channels: Six independently controllable
- Pulse shape: symmetric, bi-phasic square pulse
- Frequency: 0.3-500Hz
- Phase duration: 50-400µs, adjustable
- Amplitude: Stimulation by 'search and stimulate' probe
- Electro-acupuncture: High: 0-32mA plus/minus 25%
Low: 0-16mA plus/minus 25%
- Timer: 0-60minutes
- Output modes: (1)Constant
(2) Burst
(3) Surge
(4) Fast+Slow
(5) Sweep
(6) Random prog.1
(7) Random prog.2
(8) Random prog.3

Code
IT7001

Point Locator and Stimulator



- Four output channels
- New IC circuitry, no cross over
- Polarity reversal switch for each channel
- Sensitive point location with both light and sound indication
- Loe-Needle-Stirn treatment and location switch
- High/ Low voltage switch: low for needles, high for TENS
- Continuous, intermittent, and dense-disperse (modulated) wave forms
- Two frequency range setting: 0 to 100 and 10 to 999 Hz
- Improved alligator clip wires and detection probe
- Low battery indicator

Code

AWQ104E

PRO Acupuncture Stimulator



- Five outputs (including an output exclusively for the use of search/stimulation probe for point detection and immediate stimulation.)
- Frequency and Pulse Width are adjustable and separated into two groups, Group 1 for Ch 1 and 2, and Group 2 for Ch 3 and 4
- Two large LCD screens for displaying precise frequency of each group, and one LCD screen for showing sensitivity during acu-point detection
- Three different modes: Continuous, Modulate and Burst.
- High & Low selection of output. X1 and X5 (multiplier) of pulse rate and width adjustments
- Timer selection for different period of treatment session: Continuous, 15 minutes and 30 minutes
- Buzzer ON/OFF, for sound control during point detection
- Steady output intensity increment control, no cross over, and no interference between each output
- Music for timer alert and error warning

Code

AWQ105

Pointer Plus Locator and Stimulator



This offers sensitive acu-point or trigger point location and stimulation, which are independently adjusted with two separate control dials. The pulse frequency is fixed at 10 Hz. A flashing light and sound indicate when point detection is made.

The spring-activated detection probe allows a uniform pressure when searching. The grounding circuit is made by holding the metal of the unit or with the hand grounding pole.

- Frequency: 10Hz (fixed)
- Volts and current: 0-25 volts, 0 to 50 milliamps
- Pulse width: 220 microseconds
- Pulse shape: asymmetric biphasic square wave
- Pulse mode: continuous

Code
IT7004

Pointer Pulse



Acupuncture points locator and stimulator with infrared laser. The device measures the treatment location through sound as well as a green light indicator on the top of the unit. Location sensitivity and sound level are adjustable.

- Output channels: one
- Volts and current: Normal dial 22.5 volts, 0 to 45 milliamps Blue zone 11.0 volts, 0 to 2 microamps
- Pulse width: 260 microseconds (μ S)
- Pulse shape: asymmetric biphasic square wave with negative spike
- Pulse modes: continuous
- Laser wavelength: 650nm
- Maximum laser power: 5 mW
- Skin monitoring: sound and green light indication
- Laser operation: yellow light indication
- Power source: 9 volt alkaline battery

Code
IT7003

Pointer Pal Locator



- Easy to operate hand-held point locator. Pointer locating only, no stimulation function
- When an active point is located, a buzzing sound is produced and the pilot light on top of the unit blinks continuously

Code
IT7006

Pointer Pulse II Dual Probe Laser



- Painless and without sensation of numbness
- Harmless, no tissue damage
- Complete sterility, no-cross infection
- No heat, no adverse effect
- User-friendly design and operation

IT7003B

Pointer-Mate Acu-Locator



- The Point Mate is a pocket portable locator
- Convenient palm-size, 2.3" x 4" x 1 "
- Features sensitivity control and point location with sound and green LED light
- Includes hand ground pole and search probe with 3.5 mm jack and 9 volt battery

Code
IT7007

Moxa & Heating Devices

meridius® Manina Moxa



Meridius® Manina is a smokeless Korean moxa made from 100% real mugwort, also known as "artemesia vulgaris" or "ssuk" in Korean.

It can help reduce pain, muscle spasms and swelling and can be useful to improve blood and lymphatic circulation.

The average temperature on the skin where Meridius® Manina Moxa is applied is 43.5 degrees Celsius.

Meridius® Manina Moxa is a safe, patient-friendly product as the ash is designed to not fall off the cone. Even more, this moxa includes 40 pre-application plasters and lasts two times longer compared to standard stick-on moxa.

Code	Pieces
MA200MID	40



Stick-On Moxa

Arirang Mild Moxa

- Easy to use stick-on moxa
- Rolled cylinder of moxa on a sticky base



Code	Contents
DB203	Mild, 200pcs / box

Health Preserving Moxa

- Big sized stick-on moxa.
- Moxa rolled into a cylinder with an adhesive base



Code	Contents
DB204	25 X 25 mm, 30pcs / box

Chosei-Kyu Stick-On Moxa

- Heat regulating paper is enclosed
- Only the highest quality ibuki moxa is used



Code	Contents
IT2001W	Weak, 100pcs / box
IT2001M	Medium, 100pcs / box
IT2001S	Strong, 100pcs / box

Neo Acu-Point Moxa

- A new type of indirect moxa
- A gentle heat penetrates acupoint



Code	Contents
IT2003A	120pcs / box
IT2003B	600pcs / box

Moxa & Heating Devices

DongBang Charcoal Moxa Cone



- The unique oven baking technique produces moxa charcoal which eliminates smoke, odour associated with traditional moxa
- The cone is coated with a special compound comprised of Bio-clay, Jade, Germanium and Bio-ceramic for enhanced effect
- It requires approximately 5 minutes until effective temperature is reached
- Highly recommended to be used with DB253

Code	Contents
DB250	approx. 55pcs / container

meridius® Mini Moxa



Meridius® Mini Moxa is a traditional, all-natural compressed type of Korean moxa shaped into small tubes. It comes in two versions: smokeless and natural. The first one produces almost no smoke or odour which makes it ideal for modern-day moxa treatments, usually performed in clinics or indoor areas.

Code	Type	Contents
ME220	Natural	180 pcs / box
ME207	Smokeless	180 pcs / box

DongBang Moxa Device




- Screw type device for moxa treatment
- It has mesh in the middle of the device to hold the burning moxa, the heat passes up and down through the mesh-protected top and bottom
- The screw in the middle can adjust the height of the device controlling level of heat

Code
DB253

Smokeless Moxa

Chosei-Kyu Smokeless Moxa



- Generates no odor and little smoke
- During use, no ventilation is required, and curtains and clothes are free of odor
- A little smoke at first, but will disappear once lit
- Stable heat temperature: besides moxa, charcoal is also used as heat source

Code	Contents
IT2002S	200pcs / box
IT2002L	1000pcs / box


DongBang Smokeless Gu Gwan Moxa



- Excellent quality moxa for needle warming
- Smokeless and odor-free
- Charcoaled moxa rolled into a tube with hole in the middle
- It is to be used with the metal cap (DB206B Needle Cap).
- Insert the metal cap into the tube, place the metal cap on top of the acupuncture needle
- It can also be used with moxa devices (DB253 DONGBANG Moxa Device) and moxa burning bowls (DB209A and DB209B)
- Comes with 10 metal caps (DB206B)

Code	Contents
DB208A	200 pcs / box

Ondan Smokeless Needle Moxa



- Japanese, odour-free moxa
- Smokeless moxa
- 12mm long cylinder shape with a 3mm hole in the middle to go through the needle handle. Cylinder diameter is approx 12mm

Code	Contents
DBA2899-120	120pcs, 1 scoop & 6 baskets
DBA2899-480	480 cones

Ondan Cone Moxa



Smokeless needle moxa. Finely ground charcoal and moxa herbs are compressed into caps. The shape and form of the moxa cone is ideal for needle warming. It is smokeless, odourless and burns down to a solid ash which is easily removed from the needle head. The burning time for each piece is approximately 6 minutes under normal room temperatures. The ash remains intact even when the burning complete.

Code	Contents
DBA2901-225	225pcs
DBA2901-450	450pcs

Moxa Cap

- Needle cap to hold loose moxa
- Place the cap on the needle's handle



Code
DBA261

Basket for Ondan Moxa

- Japanese Moxa Basket for use with Ondan Moxa cones



Code	Contents
DB2899-BT-01	1pc
DB2899-BT-05	5pcs
DB2899-BT-10	10pcs

Direct Super Mini Moxa

- Pre-rolled, fine thread for direct moxa treatment
- Place the burning thread directly on the skin



Code	Contents
DB212W	0.5 mm, 450pcs / box
DB212M	0.8 mm, 350pcs / box
DB212S	1.2 mm, 300pcs / box

Metal Cap for Gu Gwan Moxa

- Insert the cap into the needle moxa
- Not suitable for Chinese-style needles



Code	Contents
DB206B	100 pcs / pack

Moxa & Heating Devices

Direct & Loose Moxa

DongBang Super Gold Moxa



The best Korean pure gold moxa

Code	Contents
DB200A	10g
DB200B	50g / bag
DB200C	100g / bag

meridius® Shinhwa Gold Moxa



Meridius® Shinhwa Mogusa Gold Moxa is a loose, premium type of Korean moxa ideally used with moxa devices. This moxa is made from 100% natural mugwort.

This herb is believed to have traditional medicinal properties that could help blood circulation and increase the general state of well-being.

Code	Contents
MEDM10	10g / box
MEDM40	40g / box

Fukuju Pure Moxa Set



Medium grade moxa for everyday professional and home use. Made from good quality ibuki moxa.

Code	Contents
IT2006	10g / box

Japanese Ibuki Gold Mountain Moxa Set



- Super pure Ibuki moxa from Japan
- Easy to mould into any shape and size for direct moxa use
- Includes 5 incense sticks

Code	Contents
IT2005	10g / box

Yomei Moxa



A needle-specific moxa, Yomei moxa features super cohesiveness for easy adhesion to the needle. As a direct moxa, it is of standard quality just below that of Fukuji Wakakusa moxa.

Code	Contents
DB2031-100	100g / box

Wakakusa Moxa



Wakakusa is the choice moxa for use in moxa needle heating (fine needle). Excellent adhesiveness with no stems. Its heat and easy lighting make it ideal for applying to acupuncture needles.

Code	Contents
IT2021-30	30g / bag
IT2021-300	300g / box

Rice Grain Shaped Direct Moxa Maker



- Very handy and practical for forming moxa rice grains
- It saves time and achieves consistent results

How to use
Place a small portion Dong Bang super pure gold moxa or ibuki moxa between the plates. Rub the plates about 10 to 20 times with light pressure back and forth. A long string of moxa is formed, from which the rice grains can be easily separated. Cork is ideal for this purpose, particularly because of its surface texture and flexibility One set of 2 boards.

Code
DB2019

Moxa Rolls

DongBang Moxa Rolls

- Longer moxa roll with super pure moxa
- Suitable for indirect moxa treatment



Code	Contents
DB205	9 x 120 mm, 20pcs / box

Kamaya Mini Smokeless Moxa

- Japanese stick-on moxa
- The adhesive on the base helps secure the moxa



Code	Contents
DB2903-120	120pcs / box
DB2903-600	600pcs / box

Hoist Smokeless Moxa Rolls

- Very condensed charcoalled moxa
- Contains a variety of herbs



Code	Contents
BN2001	5 pcs/ box

Smokeless Moxa Perforated Rolls

- Charcoalled moxa cigar with a central hole for easy lighting.



Code	Contents
BN2020	5 pcs / box

Pure Moxa Rolls



- Classical Chinese pure moxa rolls
- For indirect moxa treatment
- Long and thick rolls with pure moxa contents

Code	Contents
DBA264	10 pcs / box

Tai-Yi Moxa Rolls



- Contains mainly moxa
- Mixed with a variety of herbs (Rhizoma Et Radiz Notopterygii, Radix Aucklandiae, Radix Angelicae Dahuricae, Radix Salviae)

Code	Contents
DBA265	10 pcs / box

Moxa Roll Holder



- Moxa roll holder, holds roll safely
- Air chamber for good ventilation
- Bakelite

Code
BN2010

Gu Gwan Moxa

meridius® Gu Gwan Moxa

- Compressed moxa formed into small tube
- Features a hollow centre in the middle



Code	Contents
ME208	7 x 25 mm, 200 pcs / box

Tiger Warmer Device

- Easy to use for warming on small areas
- Warm acupoints along the meridians



Code
DB215A

Lion Warmer Device

- For powerful heat application on larger areas
- Control the heat by pressing the spring on the top



Code
DB215C

meridius® Tiger Warmer Moxa

- Compatible with Tiger Warmer Device
- Comes smokeless and natural



Code	Contents
DB215B	Natural, 5 x 80 mm, 50 pcs / box
DBA2144	Smokeless, 5 x 80 mm, 50 pcs / box

Moxa & Heating Devices

Moxa Accessories

Gu Gwan Bowl

- Ideal for indirect moxa treatment
- Features holes for acupuncture needles



Code	Type
DB209A	1 hole bowl
DB209B	5 holes bowl

OS Moxa Balm

- Can be used for direct moxa, infrared heating treatment and massage



Code	Size
MM2005	20g / jar

Shiunko Cream

- Use for direct moxa techniques, or with the tiger warmer or lion warmer



Code	Size
ITD2010S	9g / jar
ITD2010C	20g / jar

Moxa Rolls Extinguisher

- Chrome plated moxa extinguisher
- Heavy and solid



Code	Height
BN2008	40 mm

Moxa Scoop



- Ideal to remove the moxa from the needle handle
- With enclosed walls on three sides

BN2004

75 x 20 x 20 mm

Moxa Rolls Extinguisher with Plate



- Moxa extinguisher mounted in the middle of the plate
- Safe and stable to extinguish, to hold the moxa rolls and collect

Code	Diameter
BN2007	80 mm

Moxa Spoon



- Stainless steel spoon with a split in the middle
- Safe and easy to remove moxa from acupuncture needles

Code	Size
BN2011	Standard
BN2011S	Small

Moxa Torch



- Ideal to light moxa due to intense jet-flame
- Refillable with lighter gas (not fluid) available from most supermarkets
- Ergonomically designed for ease of holding
- Dispatched unfilled

Code
DBA257

Moxa Lighter



- Ideal to light moxa due to intense jet-flame
- Refillable with lighter gas (not fluid) available from most supermarkets
- Ergonomically designed for ease of holding
- Dispatched unfilled

Code
DBA258

Morning Star Japanese Incense



Japanese incense Nippon Koda in delicate natural fragrances

Code	Fragrance
DBA212AR	Amber, 50 sticks
DBA212CW	Cedarwood, 50 sticks
DBA212CM	Cinnamon, 50 sticks
DBA212GT	Green Tea, 50 sticks
DBA212JA	Jasmine, 50 sticks
DBA212LA	Lavander, 50 sticks
DBA212MU	Musk, 50 sticks
DBA212PA	Patchouli, 50 sticks
DBA212PN	Pine, 50 sticks
DBA212RO	Rose, 50 sticks
DBA212SW	Sandalwood, 50 sticks
DBA212VN	Vanilla, 50 sticks
DBA212YZ	Yuzu, 50 sticks
DBA212IR	Iris, 50 sticks
DBA212LO	Lotus, 50 sticks
DBA212FG	Fig, 50 sticks
DBA212GA	Gardenia, 50 sticks
DBA212MM	Mimosa, 50 sticks



Maeil Koran Incense



Sandlewood scented incence

Code	Type	Contents
DB210A	Smokeless	250 pcs / box
DB210B	Natural	200 pcs / box

Moxa & Heating Devices

HerbaChaud



- Enhanced by a traditional Chinese blend of herbs
- Promotes general physical and mental wellbeing
- Has a pleasant and relaxing effect on the body
- Soothes pain and stops inflammation
- Used for the local relief of muscular, rheumatic and arthritic pain, including tense muscles (back and neck), contusions, strains, and tennis elbow

Code	Contents
CT2001S	2 pcs / pack
CT2001L	6 pcs / pack

PREMIO 10 Moxa



- Has an infrared emitter which is equivalent to emission of burning moxa
- This Premio moxa is innovative, safe, effective and pleasant

Code
AC2001

Boisun Ear Candles



- Original "Biosun" Hopi Ear Candles. Treatment for ear problems such as: aches, infections and wax build up

Code	Contents
HP515-01	1 pair / pack
HP515-05	5 pairs / pack

Heating Lamps

TDP Infrared Heating Lamps



TDP is the abbreviation of 'the special electromagnetic spectrum' in Chinese. It has been frequently used in clinics and is also referred to as the 'Miracle Lamp'.

The TDP lamp has a plate coated with 33 important minerals for the human body. This charged thermal plate gives a spectrum between 2-50 Micro in wavelength, and between 28-34 mw/ cm2 intensity. It also gives heat element resulting in a special form of electromagnetic waves, in which both wavelength and strength are ideal for the human body.

As a result, the lamp works like moxa, giving a heating affect to depth of 3-5 cm and it is found to be very effective for patients. It can be used immediately following acupuncture and without needles, giving a similar result to the use of warm needles.

Code	Contents
HT2012UK	Table-top with digital timer
HT2000-165	Replacement plate, 16.5cm
HT2000-125	Replacement plate, 12.5cm
HT2036UK	Single Head TDP Lamp
HT2033UK	Double Head TDP Lamp

Meridian & Muscle Model

- Half of the model displays meridian and points, the other half shows the muscle



Code	Size
DBA6021	80 cm
DBA6022	65 cm

Meridian Model

- Soft white vinyl with 360 points and 14 meridians
- Both Chinese and English coded points illustration



Code	Type
DBA601	Male, 50 cm
DBA602	Female, 48 cm

Hand Model

- Plastic model shows number coded acupuncture points with a booklet



Code	Size
DBA605A	22 cm
DBA605B	13 cm

Foot Model

- Made of soft vinyl
- Illustrates the acupuncture points



Code	Size
DBA607A	15 cm

Foot Model Pair of Left and Right



- A pair of soft vinyl feet that shows reflex zones of the foot
- Ideal for reflexology

Code	Size
DBA605F2	18 x 8 cm

Ear Model



- Made of soft vinyl
- Rests on a plastic base
- Both Chinese and English coded illustration
- Coding booklet include

Code	Size
DBA604A	22 cm
DBA6048	13 cm

Head Model



- Made of soft vinyl
- Features a plastic base
- Both Chinese and English coded illustration
- With a booklet
- Height 22cm

Code
DBA605

Silicone Realistic Ear Model



- Life-size high quality silicone ear models
- Includes pair of two: one left ear and one right ear
- No visible holes form after practising with needles

Code	Type
DAV1000	Pair of 2

Animal Models



- Plastic animal model mounted on a plastic stand
- Front side of model shows acupuncture points numbered
- Comes with a booklet describing point locations and recommendations for treatments

Code	Size	Type
DBA608	26 x 7 x 23 cm	Horse
DBA609	30 x 9 x 28 cm	Dog
DBA610	28 x 8 x 19 cm	Cat
DBA612	26 x 8 x 14 cm	Cow

Dr Life Intermittent Pneumatic Compression Machines



These Intermittent Pneumatic Compression (IPC) machines are made by the world-renowned South Korean firm, Daesung Maref. Established in 1986, Daesung Maref now provide 90% of the IPC sales in the South Korean market and export to 30 countries worldwide.

With Intermittent Pneumatic Compression, the limb is surrounded by an inflatable cuff with several chambers, which are inflated sequentially so that there is a wave of compression travelling up the limb in short bursts. Venous emptying is made easy and an active circulation is expected. It is also useful in the arm to reduce lymphoedema after mastectomy.

Benefits

- Compact and simple design for home use
- Restful and calming
- Helps relieve stress for elderly and hospitalized patients
- Waist cuff can help to restore normal bowel functioning
- Helps to prevent DVT, varicose veins, lymphoedema, numbness of the extremities, poor venous return

References

- Arcelus JI et al, Home use of impulse compression of the foot and compression stocking in the treatment of Chronic venous insufficiency J Vase Surg 2007 ;34:805-11
- Kalodiki E et al, Intermittant Pneumatic Compression (IPC) in the Treatment of Peripheral Arterial Occlusive Disease (PAOD) - A Useful Tool or Just Another Device Eur J Vase Endovasc Surg 2007; 33, 309-37
- van Bemmelen PS et al, rapid Intermittent compression increases skin circulation in chronically ischaemic legs with infra-popliteal arterial obstruction. Vasa 2000; 29: 47-52
- Labropoulos N et al, Haemodynamic effects of intermittent pneumatic compression in patients with critical limb ischaemia. J Vase Surg 2005; 42: 710 - 716

Code	Type
LX7	Dr. Life Machine
LX9	Dr. Life Machine
DL2003	Extention zipper



Luxury Embossed Couch Rolls



- 9 rolls per case
- Full 20"/508mm wide roll
- Premium quality 2ply tissue
- Easy tear perfect perforations
- Bright white in colour
- 32m per roll
- Embossed for softness and strength
- 100% pure paper
- High absorbency
- Enhanced texture
- Increased client comfort
- Made in the UK

Code
DBA6058809

Pre Injection Swabs



- Non-woven fabric swab saturated with 70% isopropyl alcohol for skin cleansing prior to administering an injection
- Sterile and individually wrapped
- 100 pcs per box

Code	Contents
DBA605108S	1 box
DBA605108S-10	10 boxes

Mini Pre-Surgical Marking Pen



- Mini Pre-Surgery Skin Marker with Scrub Resistant Ink
- For best results, degrease skin before use with alcohol wipe
- Mark skin, allow to dry

Code
DBA605100

Hard Surface Wipes



- 70% Isopropyl alcohol solution for fast bactericidal action
- Disinfects surfaces quickly and effectively
- Effective against MRSA, E-Coli and a range of micro organisms

Code	Contents
DBA6051077	200 pcs / tube

Treatment Room Equipment

Power Free Latex Gloves



Code	Size
DBA6051240	Extra Small
DBA6051210	Small
DBA6051220	Medium
DBA6051230	Large

Sharps Disposal Box



- Material: Polypropylene
- Ideal for travelling, easily fitted into a bag or glove compartment
- Non-return valve

Code	Size
DBA60510-05	0.5 litre
DBA60510-10	1 litre
DBA60510-25	2.5 litre
DBA60510-375	3.75 litres
DBA60510-5	5 litres



Orders may placed online, email or by phone.

Online: www.eu.meridiusmedical.com

Email: eu@meridiusmedical.com

Phone: +49 (0) 208 20580580

Address: Heiderhöfen 23, Oberhausen, 46049, Germany

VAT no:

Terms & Conditions

All products shown are subject to availabilty. We reserve the right to refuse to supply to any individual or company.

E-catalogue can be downloaded from www.eu.meridiusmedical.com.

Meridius Medical GmbH

Heiderhöfen 23
Oberhausen,
46049, Germany

www.eu.meridiusmedical.com

Email: eu@meridiusmedical.com

Tel: +49 (0) 208 20580580